

Table 1. Primers and probes used for Real-time PCR

Gene	Primer and probe sequence
<i>GAPDH</i>	Forward 5- ACACCCACTCCTCCACCTTTG-3 Reverse 5- GTCCACCACCCTGTTGCTGTA-3
<i>p35</i>	Forward 5-AACCAAGAATGAGAGTTGCCTAA- Reverse 5-GAACTCCACCTGGTACATCTTC-3
<i>IL-37</i>	Forward 5-AGTGAGGTCAGCGATTAGGAA-3 Reverse 5-TTTAGTGAGCAGGTTTGGTGTTC-3
<i>EBI</i>	Forward 5-GCCTTTCATAACAGAGCACATCAT-3 Reverse 5-AACGGATCCAGTACTTCAGTGAG-3