

Table S1. Primers used for Real-Time polymerase chain reaction assay

Gene	Primers (5'-3')	Product Size (bp)
<i>Galectin-9</i>	F: TTTCTGGGACTATTCAAGGAG R: GAAGTGGAAGGCAATGTCA	137 bp
<i>PD-L1</i>	F: CTATGGTGGTGCCGACTACAA R: CTGCTTGTCAGATGACTTCG	159 bp
<i>CD155</i>	F: GGACGGCAAGAATGTGAC R: CCAGTTGTTATCATAGCCAGAG	124 bp
<i>β-actin</i>	F: CCTTCCTGGGCATGGAGTCCT R: TGGGTGCCAGGGCAGTGAT	174 bp

F, Forward primer; R, reverse primer